

SINCE 1908  
**wessels**  
 company

**SUBMITTAL**

**TYPE:    DUPLEX GMP SERIES  
 GLYCOL MAKE-UP PACKAGE**


Submittal Sheet No. M-3223      Date 7/01

JOB _____  Unit Tag No. _____ Engineer _____ Contractor _____	Wessels Representative _____ _____ _____ Order No. _____      Date _____ Submitted By _____      Date _____ Approved By _____      Date _____
---	--

**DUPLEX GLYCOL MAKE-UP PACKAGE**

Fully automated and autonomous. Complete with: low level cut-off and alarm arrangement including a 110V signal for remote alarm, isolation valves, strainer, pressure tank with pressure control, pressure reducing valve and gauge, translucent polyethylene solution container with lid designed to accommodate relief valve piping, magnetic starter, 110V, 60 Hz motor and controls. All mounted on a sturdy steel frame with 1/2" system connection and factory preset 12 psig (field adjustable) discharge.

**GMPD-23050 & 23100 - 1/3 HP**  
**GMPD-25050 & 25100 - 1/2 HP**



Model Number	Make-Up Capacity (gpm @ psi)	Pressure Range (psi)	Solution Container (gal)	Dimensions (in inches)		Approx. Weight (lbs)
				A	B	
GMPD-23050	1.8 @ 50	10-50	50	52	34	200
GMPD-23100	1.8 @ 50	10-50	100	76	34	216
GMPD-25050	1.7 @ 70	10-70	50	52	34	206
GMPD-25100	1.7 @ 70	10-70	100	76	34	222

NOTE: Please consult factory for other sizes and/or models.

QUANTITY      MODEL  
 \_\_\_\_\_ GMPD- \_\_\_\_\_